

Il contenuto non è disponibile nella lingua selezionata.

Verrà mostrato il contenuto in un'altra delle lingue disponibili. Il tuo browser potrebbe avere delle funzionalità di traduzione.

GPIB Instrument Control Tutorial

Updated Apr 17, 2023

Environment

Hardware

- [GPIB Instrument Control Device](#)

This tutorial shows you [how to set up National Instrument hardware and software to communicate with a GPIB instrument](#) in LabVIEW.

- Learn more about [GPIB communication](#)
- Learn more about [Software Required for Instrument Control \(GPIB, Serial, VISA, USB, etc\)](#)

Prerequisites

Consult the pinout or documentation for your GPIB instrument and [make sure you have the correct cables and GPIB instrument control device to connect the instrument to the computer.](#)

Cables

- [Shop NI GPIB cables](#)
- [GPIB Cabling Requirements and Limitations](#)

GPIB Instrument Control Device

- [Shop GPIB instrument control devices](#)

Install NI Software and Drivers

[Install the latest compatible version of LabVIEW](#), the [NI-VISA](#) driver, and the [NI-488.2](#) driver. Refer to the following for compatibility information.

- [LabVIEW and Microsoft Windows Compatibility](#)
- [NI-VISA and Operating System Compatibility](#)
- [NI-VISA and LabVIEW Version Compatibility](#)
- [NI-488.2 Operating System Support](#)
- [NI-488.2 Supported Versions for GPIB Devices and Modules](#)
- [NI-488.2 and LabVIEW Version Compatibility](#)

Install GPIB Controller

If you are using an NI GPIB controller or a NI GPIB-ENET device, refer to the [GPIB Hardware Installation Guide](#) to install the hardware.

Configure Hardware in MAX

In Measurement and Automation Explorer (MAX), NI-VISA detects the instrument and creates the VISA resource that you use to communicate with the instrument.

1. Open MAX from the Start menu.
2. Expand **Devices and Interfaces**. You should see the VISA resource name of the NI GPIB controller listed as GPIB<controller>::INSTR where <controller> is the controller number. If not specified, the controller number defaults to 0.

If you are unable to see the controller, refer to [Cannot Detect NI GPIB-USB HS Device in MAX](#).

3. After the controller has been successfully recognized, securely connect the GPIB instrument to the controller via a GPIB cable and power it on.
4. In MAX, select the GPIB controller and select **Scan for Instruments**. NI-VISA will attempt to identify the connected GPIB instrument by conducting a FindLstn() test and a *IDN? query. The instrument VISA resource name should appear under the GPIB controller resource as GPIB::<primary address>::INSTR.

Note: If your instrument is not SCPI compliant, it will not appear in MAX. Check the instrument's user manual for supported commands and continue to the Test Communication section below.

If **Scan for Instruments** is missing, refer to [NI MAX GPIB Scan for Instruments Missing](#).

If you experience issues with **Scan for Instruments**, refer to [Can't See GPIB Instrument Using Scan for Instruments in MAX](#).




Figure 1. Scanning for Instruments in MAX

5. After your instrument has been successfully recognized, select the instrument in MAX and then the VISA Properties tab at the bottom of the window. You can configure the VISA Alias of the instrument for easier identification.

Test Communication

After set-up and configuration, use an interactive control to verify communication with the instrument and test different instrument commands.

If your device is SCPI compliant, select its VISA resource and select **Communicate with Instrument**.

1. The NI-488.2 Communicator should launch.
2. Type a supported command in the **Send String** field. Refer to manufacturer documentation for a list of valid commands.
3. Select one of the following:

Query – writes the command to the instrument and reads the response back.

Write – writes the command to the instrument.

Read – reads the response from the instrument.

If your device is not SCPI compliant, use GPIB Interactive Control (IBIC) for more advanced communication.

1. In MAX, navigate to the **Tools** tab >> **NI-488.2 >> Interactive Control**.
2. Enter help to launch IBIC help documentation and refer to your instrument's manual for supported commands.

If you are able to communicate with the instrument in IBIC, then you should also be able to communicate with the instrument from any other programming language, such as G code or C.

Install Instrument Driver

After completing your configuration and verifying communication with the instrument, you need an easy way to transition from interactive mode to programming mode. The best way to achieve this is by using instrument drivers.

Refer to [Choosing the Right Interface to Control Instruments in LabVIEW](#) to determine which instrument driver type is best for your application. Follow the links in the document to install and use the instrument driver.

Run an Example

The instrument driver, NI-488.2, and NI-VISA driver will install examples and the API required to communicate with your instrument. Use these examples as programming references and an easy way to start using your instrument. For more information on whether to use the NI-488.2 or NI-VISA API, refer to [NI-VISA Help: GPIB Summary](#).

Note: Close MAX before you start running a LabVIEW example. Having MAX Test Panels open and running a LabVIEW example simultaneously will result in an access resource error.

Instrument Driver Examples

If you installed an instrument driver, refer to the following articles for instructions on running examples:

- [Using Plug-and-Play Instrument Drivers in LabVIEW](#)
- [Using IVI Drivers in LabVIEW or LabWindows™/CVI™](#)

NI-488.2 or NI-VISA Driver Example

If you are using NI-488.2 or NI-VISA, refer to the shipping examples Simple GPIB.vi. or GPIB with VISA Functions.vi.

1. To locate the examples, open LabVIEW and navigate to the **Help** tab >> **Find Examples >> Hardware Input and Output >> GPIB**.
2. Ensure the proper VISA resource name is selected.
3. Run the VI.

For more information on the NI-488.2 functions, refer to [LabVIEW Help: GPIB Functions](#).

For more information on the NI-VISA functions, refer to [LabVIEW Help: VISA Functions and VIs](#).

For a comparison on the NI-488.2 and NI-VISA functions, refer to [NI-VISA Help: Comparison Between NI-VISA and NI-488.2 APIs](#).

Related Links

More GPIB Resources

- [NI-488.2 User Manual](#)
- [NI VISA Programmer Reference Manual](#)
- [Using Keysight \(HP Agilent\) GPIB Devices with NI Software](#)
- [Using Both NI GPIB and Keysight HPIB Controllers in the Same System](#)
- [Computer to Computer Communication Via GPIB](#)

- [NI GPIB Cabling Requirements and Limitations](#)

Other Instrument Control Tutorials

- [Ethernet Instrument Control Tutorial](#)
- [Serial Instrument Control Tutorial](#)
- [USB Instrument Control Tutorial](#)

Other Support Options

Ask the NI Community

Collaborate with other users in our discussion forums

- [Search the NI Community for a solution](#)

Request Support from an Engineer

A valid service agreement may be required, and support options vary by country.

- [Open a service request](#)
- [Purchase or renew support services](#)